


Anti-Pak6 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nqu5

Item Number: CM0030803



Introduction

Anti-PAK6 rabbit polyclonal antibody validated for WB and ELISA. Detects human and mouse PAK6 (also known as PAK5), a serine/threonine kinase involved in transcription regulation. Selectively expressed in brain and testis. Ideal for research use.

[Learn More](#)

Parameter	Value
Product Name	PAK6 Rabbit pAb
Remarks/Alias	PAK5; PAK6
Species	Human
GeneID (Human)	56924
GeneID	56924
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human PAK6 (NP_064553.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NQ5
Protein Weight	75 kDa
Shipping	Ice bag