

Anti-Palb2 Rabbit Polyclonal Antibody For Wb And Elisa - Q86Yc2

Item Number: CM0030701



Introduction

PALB2 Rabbit Polyclonal Antibody validated for WB and ELISA. This antibody targets human Partner and localizer of BRCA2, a key protein in homologous recombination repair and Fanconi anemia pathway. Ideal for cancer research.

[Learn More](#)

Parameter	Specification
Product Name	PALB2 Rabbit pAb
Remarks/Alias	FANCN; PNCA3; PALB2
Species	Human
GeneID (Human)	79728
GeneID	79728
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 99-198 of human PALB2 (NP_078951.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q86YC2
Protein Weight	131kDa
Shipping	Ice bag