

Anti-Pald1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Ule6

Item Number: CM0017865

Introduction

Rabbit polyclonal antibody against human PALD1 (paladin), raised against an N-terminal recombinant protein. Validated for WB, IF/ICC, and ELISA in human, mouse, and rat samples; detects the ~97 kDa cytosolic protein.

[Learn More](#)

Parameter	Value
Product Name	PALD1 Rabbit pAb
Remarks/Alias	PALD; KIAA1274; PALD1
Species	Human
GeneID (Human)	27143
GeneID	27143
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human PALD1 (NP_055246.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9ULE6
Protein Weight	97kDa
Shipping	Ice bag