

Anti-Pam16 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y3D7

Item Number: CM0030877

Introduction

Rabbit polyclonal antibody against human PAM16 (TIM16, MAGMAS), validated for Western blot and ELISA; cross-reacts with mouse. Detects a 14 kDa protein; suitable for mitochondrial import research.

[Learn More](#)

Parameter	Value
Product Name	PAM16 Rabbit pAb
Remarks/Alias	TIM16; MAGMAS; SMDMDM; TIMM16; CGI-136; PAM16
Species	Human
Gene ID (Human)	51025
Gene ID	51025
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-125 of human PAM16 (NP_057153.8)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y3D7
Protein Weight	14kDa
Shipping	Ice bag