


Anti-Pan 14-3-3 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P27348 / P31946 / P61981 / P62258 / P63104 / Q04917 / P31947

Item Number: CM0006354



Introduction

High-quality rabbit monoclonal antibody recognizing all seven human 14-3-3 isoforms (θ , β , γ , ϵ , ζ , η , σ). Validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat samples.

[Learn More](#)

Parameter	Value
Product Name	pan 14-3-3 Rabbit mAb
Remarks/Alias	1C5; HS1; 14-3-3; GW128; YWHAA; KCIP-1; HEL-S-1; HEL4; YWHAD; HEL-S-3; POPCHAS; HEL-S-93; 14-3-3-zeta; pan 14-3-3
Species	Human
GenelD (Human)	10971/7529/7532/7531/7534/7533/2810
GenelD	10971/7529/7532/7531/7534/7533/2810
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids human 1-100 of human pan 14-3-3 (NP_006817.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P27348/P31946/P61981/P62258/P63104/Q04917/P31947
Protein Weight	27kDa/28KDa/29KDa
Shipping	Ice bag