


Anti-Pan 14-3-3 Rabbit Polyclonal Antibody For Wb, Elisa - P27348 / P31946 / P61981 / P62258 / P63104 / Q04917 / P31947

Item Number: CM0022017



Introduction

Rabbit polyclonal antibody targeting pan 14-3-3 (YWHAB, YWHAZ, etc.), validated for WB and ELISA. Cross-reacts with human, mouse, rat, and African green monkey. Ideal for cell signaling and neurodegeneration research.

[Learn More](#)

Parameter	Value
Product Name	pan 14-3-3 Rabbit pAb
Remarks/Alias	YWHA; S1C5; HS1; 14-3-3; pan 14-3-3
Species	Human
GeneID (human)	10971/7529/7532/7531/7534/7533/2810
GeneID	10971/7529/7532/7531/7534/7533/2810
Immunogen	Recombinant protein of human pan 14-3-3
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat, African green monkey
SWISS	P27348/P31946/P61981/P62258/P63104/Q04917/P31947
Protein Weight	27kDa/28kDa
Shipping	Ice bag