

Anti-Pan Cadherin Polyclonal Antibody For Wb, Ip, Elisa - P12830 / P19022 / P22223 / P55283

Item Number: CM0022697



Introduction

CamelBio Pan Cadherin Rabbit Polyclonal Antibody for WB, IP, ELISA in human, mouse, rat. Recognizes conserved cadherin region for broad cell adhesion and diagnostic research use.

[Learn More](#)

Parameter	Value
Product Name	Pan Cadherin Rabbit pAb
Species	Human
GeneID (Human)	999/1000/1001/1002
GeneID	999/1000/1001/1002
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 800-882 of human pan-cadherin (NP_004351.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12830/P19022/P22223/P55283
Shipping	Ice bag