

## Anti-Pan-Trk Polyclonal Antibody For Wb And Elisa - P04629/Q16620/Q16288

Item Number: CM0021455



### Introduction

High-quality rabbit polyclonal antibody targeting human pan-Trk (NTRK1/2/3), validated for Western blot (WB) and ELISA with mouse cross-reactivity. Ideal for neurotrophin signaling and cancer research.

[Learn More](#)

| Parameter        | Value                                                                                                   |
|------------------|---------------------------------------------------------------------------------------------------------|
| Product Name     | pan-Trk Rabbit pAb                                                                                      |
| Species          | Human                                                                                                   |
| Gene ID (Human)  | 4914/4915/4916                                                                                          |
| Gene ID          | 4914/4915/4916                                                                                          |
| Immunogen        | Synthetic peptide corresponding to a sequence within amino acids 697-796 of human pan-Trk (NP_002520.2) |
| Source           | Rabbit                                                                                                  |
| Category         | Polyclonal Antibodies                                                                                   |
| Application      | WB, ELISA                                                                                               |
| Cross Reactivity | Mouse                                                                                                   |
| SWISS            | P04629/Q16620/Q16288                                                                                    |
| Shipping         | Ice bag                                                                                                 |