


Anti-Pan-Trk Rabbit Monoclonal Antibody For Wb, Elisa - P04629 / Q16620 / Q16288

Item Number: CM0007379



Introduction

Pan-Trk rabbit monoclonal antibody validated for WB and ELISA applications. Detects human, mouse, rat TrkA/B/C receptors. Ideal for neurotrophin signaling research and neuronal differentiation studies.

[Learn More](#)

Parameter	Value
Product Name	pan-Trk Rabbit mAb
Species	Human
Gene ID	4914, 4915, 4916
Immunogen	A synthesized peptide derived from human pan-Trk
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04629, Q16620, Q16288
Shipping	Ice bag