


Anti-Pannexin-1 Polyclonal Antibody For Wb, Elisa - Q96Rd7

Item Number: CM0027498



Introduction

High-quality rabbit polyclonal antibody targeting human Pannexin-1 (PANX1), suitable for WB and ELISA. Cross-reacts with mouse. Validated for research on ion channels, apoptosis, and calcium signaling. Swiss-Prot: Q96RD7.

[Learn More](#)

Parameter	Value
Product Name	Pannexin 1 Rabbit pAb
Remarks/Alias	PX1; MRS1; OOMD7; OZEMA7; UNQ2529; Pannexin 1
Species	Human
Gene ID (Human)	24145
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-150 of human Pannexin 1 (NP_056183.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96RD7
Protein Weight	48kDa
Shipping	Ice bag