

Anti-Pap Rabbit Polyclonal Antibody For Wb, Elisa - P51003

Item Number: CM0026312



Introduction

High-quality rabbit polyclonal antibody targeting human poly(A) polymerase alpha (PAPOLA/PAP). Validated for Western blot and ELISA, suitable for studying mRNA polyadenylation. Ideal for diagnostic and research applications.

[Learn More](#)

Parameter	Value
Product Name	PAP Rabbit pAb
Remarks/Alias	PAP; PAP-alpha
Species	Human
GeneID (Human)	10914
GeneID	10914
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 605-745 of human PAP (NP_116021.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P51003
Protein Weight	83kDa
Shipping	Ice bag