

Anti-Papd4 (Gld2/Tut2) Polyclonal Antibody For Wb And Elisa - Q6Piy7

Item Number: CM0015800

Introduction

Rabbit polyclonal antibody targeting human PAPD4 (GLD2/TUT2), validated for WB and ELISA in human, mouse, and rat samples. Suitable for studying cytoplasmic polyadenylation and miRNA regulation.

[Learn More](#)

Parameter	Value
Product Name	PAPD4 Rabbit pAb
Remarks/Alias	APD4; GLD2; TUT2; PAPD4
Species	Human
GenelD (human)	167153
GenelD	167153
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human PAPD4 (NP_776158.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6PIY7
Protein Weight	56kDa
Shipping	Ice bag