

## Anti-Papd7 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q5Xg87

Item Number: CM0015724

### Introduction

Rabbit polyclonal antibody against human PAPD7/TENT4A, validated for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Target involved in mRNA stabilization and RNA quality control.

[Learn More](#)

| Parameter        | Value                                                                                                                                   |
|------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | PAPD7 Rabbit pAb                                                                                                                        |
| Remarks/Alias    | LAK1; POLK; POLS; TRF4; LAK-1; PAPD7; TRF41; TRF4-1; TUTASE5                                                                            |
| Species          | Human                                                                                                                                   |
| GeneID (Human)   | 11044                                                                                                                                   |
| GeneID           | 11044                                                                                                                                   |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 303-542 of human PAPD7 (NP_008930.1). |
| Source           | Rabbit                                                                                                                                  |
| Category         | Polyclonal Antibodies                                                                                                                   |
| Application      | WB, IHC-P, ELISA                                                                                                                        |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                       |
| SWISS            | Q5XG87                                                                                                                                  |
| Protein Weight   | 85kDa                                                                                                                                   |
| Shipping         | Ice bag                                                                                                                                 |