


Anti-Papss1 Polyclonal Antibody For Wb - 043252

Item Number: CM0015698



Introduction

Rabbit polyclonal antibody to human PAPSS1 (SK1/ATPSK1). Validated for WB and ELISA, reacting with human, mouse, rat. Ideal for sulfate activation pathway and sulfotransferase research.

[Learn More](#)

Parameter	Value
Product Name	PAPSS1 Rabbit pAb
Remarks/Alias	SK1; PAPSS; ATPSK1; PAPSS1
Species	Human
GeneID (human)	9061
GeneID	9061
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 190-260 of human PAPSS1 (NP_005434.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43252
Protein Weight	71kDa
Shipping	Ice bag