

Anti-Par2 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P55085

Item Number: CM0010308



Introduction

Recombinant rabbit monoclonal antibody against human PAR2 (P55085), validated for WB, IHC-P, ELISA. Suitable for studies of inflammation, immunity, and endothelial barrier function.

[Learn More](#)

Parameter	Value
Product Name	PAR2 Rabbit mAb
Remarks/Alias	PAR2; GPR11
Species	Human
Gene ID	2150 (Human)
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 298-397 of human PAR2 (P55085)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55085
Protein Weight	44kDa
Shipping	Ice bag