

Anti-Par2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P55085

Item Number: CM0025996



Introduction

Validate PAR2 (F2RL1/GPR11) expression in human, mouse, rat with this rabbit polyclonal antibody. Ideal for WB, IF/ICC, ELISA. Explore inflammatory and barrier signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	PAR2 Rabbit pAb
Remarks/Alias	PAR2; GPR11
Species	Human
GeneID (Human)	2150
GeneID	2150
Immunogen	Recombinant fusion protein containing a sequence within amino acids 15-94 of human PAR2 (NP_005233.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55085
Protein Weight	44kDa
Shipping	Ice bag