

Anti-Par4 Rabbit Polyclonal Antibody For Wb, Elisa - Q96Iz0

Item Number: CM0027204



Introduction

High-quality anti-PAR4 rabbit polyclonal antibody for Western blot and ELISA, detecting human, mouse, and rat PAR4. Ideal for apoptosis and cancer research. - Q96IZ0

[Learn More](#)

Parameter	Value
Product Name	PAR4 Rabbit pAb
Remarks/Alias	PAR4; Par-4
Species	Human
Gene ID (Human)	5074
Gene ID	5074
Immunogen	Synthetic peptide; A synthetic peptide corresponding to a sequence within amino acids 200-300 of human PAR4 (NP_002574.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96IZ0
Protein Weight	37kDa
Shipping	Ice bag