


Anti-Parathyroid Hormone (Pth) Monoclonal Antibody For Wb, If-P, Elisa - P01270

Item Number: CM0005538



Introduction

Rabbit monoclonal antibody targeting human Parathyroid Hormone (PTH), suitable for WB, IF-P, and ELISA. Cross-reacts with mouse and rat. Ideal for calcium homeostasis and bone metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	Parathyroid Hormone (PTH) Rabbit mAb
Remarks/Alias	FIH1; PTH1; Parathyroid Hormone (PTH)
Species	Human
Gene ID (Human)	5741
Gene ID	5741
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Parathyroid Hormone (PTH) (P01270)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01270
Protein Weight	13kDa
Shipping	Ice bag