

Anti-Parathyroid Hormone (Pth) Rabbit Polyclonal Antibody For Wb, Elisa - P01270

Item Number: CM0026761



Introduction

This rabbit polyclonal antibody targets human PTH (SWISS P01270) and cross-reacts with mouse and rat. Suitable for WB and ELISA, it aids in calcium metabolism research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	Parathyroid Hormone (PTH) Rabbit pAb
Remarks/Alias	FIH1; PTH1; Parathyroid Hormone (PTH)
Species	Human
Gene ID (Human)	5741
Gene ID	5741
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-115 of human Parathyroid Hormone (PTH) (NP_000306.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01270
Protein Weight	13kDa
Shipping	Ice bag