

Anti-Pard6A Rabbit Polyclonal Antibody For Wb And Elisa - Q9Npb6

Item Number: CM0027272



Introduction

High-quality anti-PARD6A rabbit polyclonal antibody for WB and ELISA. Validated for human and mouse samples, targeting the PAR6C adapter protein involved in cell polarity and tight junctions.

[Learn More](#)

Parameter	Value
Product Name	PARD6A Rabbit pAb
Remarks/Alias	PAR6; PAR6C; TAX40; PAR-6A; TIP-40; PAR6alpha; PARD6A
Species	Human
Gene ID (Human)	50855
Gene ID	50855
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human PARD6A (Q9NPB6)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NPB6
Protein Weight	37kDa
Shipping	Ice bag