

Anti-Parg Rabbit Monoclonal Antibody For Wb, Elisa - Q86W56

Item Number: CM0004611



Introduction

Anti-PARG rabbit monoclonal antibody validated for WB and ELISA; specifically targets human poly(ADP-ribose) glycohydrolase. Ideal for studying ADP-ribosylation, DNA damage response, and chromatin remodeling pathways.

[Learn More](#)

Parameter	Value
Product Name	PARG Rabbit mAb
Remarks/Alias	PARG99; PARG
Species	Human
Gene ID (Human)	8505
Immunogen	Recombinant protein corresponding to a sequence within amino acids 777-976 of human PARG (NP_003622.2)
Source	Rabbit
Category	Monoclonal Antibody
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q86W56
Protein Weight	111kDa
Shipping	Ice bag