

Anti-Parn Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O95453

Item Number: CM0024851



Introduction

Rabbit polyclonal antibody against human PARN, a poly(A)-specific ribonuclease. Validated for WB, IF/ICC, ELISA. Recombinant immunogen (aa 1-280). Ideal for mRNA decay and post-transcriptional regulation studies.

[Learn More](#)

Parameter	Value
Product Name	PARN Rabbit pAb
Remarks/Alias	DAN; DKCB6; PFBMFT4; PARN
Species	Human
GeneID (Human)	5073
GeneID	5073
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human PARN (NP_002573.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	O95453
Protein Weight	73kDa
Shipping	Ice bag