


Anti-Parp1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P09874

Item Number: CM0012831



Introduction

The PARP1 Rabbit Polyclonal Antibody detects human, mouse, and rat PARP1 (113kDa) in WB, IHC-P, IF/ICC, and ELISA. It targets the C-terminal region (aa 662-1014) and is ideal for DNA damage repair, apoptosis, and transcription studies.

[Learn More](#)

Parameter	Value
Product Name	PARP1 Rabbit pAb
Remarks/Alias	PARP; PARS; PPOL; ADPRT; ARTD1; ADPRT1; PARP-1; ADPRT 1; pADPRT-1; Poly-PARP; PARP1
Species	Human
GeneID (Human)	142
GeneID	142
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 662-1014 of PARP1 (NP_001609.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09874
Protein Weight	113kDa
Shipping	Ice bag