


Anti-Parp1 Rabbit Monoclonal Antibody For If/Icc, If-P, Ihc-P, Elisa - P09874

Item Number: CM0003758



Introduction

Anti-PARP1 rabbit monoclonal antibody validated for IF/ICC, IF-P, IHC-P, ELISA. Cross-reacts with human, mouse, and rat PARP1. Essential for DNA repair, apoptosis, and transcription research.

[Learn More](#)

Parameter	Value
Product Name	PARP1 Rabbit mAb
Remarks/Alias	PARP; PARS; PPOL; ADPRT; ARTD1; ADPRT1; PARP-1; ADPRT 1; pADPRT-1; Poly-PARP; PARP1
Species	Human
Gene ID (Human)	142
Gene ID	142
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 700-800 of human PARP1 (P09874).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09874
Protein Weight	113kDa
Shipping	Ice bag