


Anti-Parp1 Rabbit Monoclonal Antibody For Wb, If, Ihc, Elisa - P09874

Item Number: CM0007772



Introduction

Rabbit monoclonal antibody targeting human PARP1. Validated for Western blot, immunofluorescence, immunohistochemistry, and ELISA. Reacts with human, mouse, and rat. Suitable for DNA repair, apoptosis, and transcription research.

[Learn More](#)

Parameter	Specification
Product Name	PARP1 Rabbit mAb
Remarks/Alias	PARP; PARS; PPOL; ADPRT; ARTD1; ADPRT1; PARP-1; ADPRT 1; pADPRT-1; Poly-PARP
Species	Human
Gene ID (Human)	142
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 200-300 of human PARP1 (NP_001609.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09874
Protein Weight	113kDa
Shipping	Ice bag