

Anti-Parp12 Polyclonal Antibody For Wb, Elisa - Q9H0J9

Item Number: CM0018108



Introduction

Rabbit polyclonal antibody targeting human PARP12, validated for WB and ELISA. Reacts with human, mouse, and rat. Ideal for research on antiviral responses, Golgi transport, and mitochondrial function.

[Learn More](#)

Parameter	Value
Product Name	PARP12 Rabbit pAb
Remarks/Alias	ZC3H1; ARTD12; MST109; MSTP109; ZC3HDC1; PARP12
Species	Human
Gene ID (Human)	64761
Gene ID	64761
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 462-701 of human PARP12 (NP_073587.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H0J9
Protein Weight	79kDa
Shipping	Ice bag