

Anti-Parp2 Rabbit Monoclonal Antibody For Wb And Elisa - Q9Ugn5

Item Number: CM0004331



Introduction

Recombinant anti-PARP2 rabbit monoclonal antibody validated for western blot and ELISA, targeting human PARP2 (Q9UGN5) for DNA repair and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	PARP2 Rabbit mAb
Remarks/Alias	ARTD2; ADPRT2; PARP-2; ADPRTL2; ADPRTL3; pADPRT-2; PARP2
Species	Human
Gene ID (Human)	10038
Gene ID	10038
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 484-583 of human PARP2 (Q9UGN5)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9UGN5
Protein Weight	66kDa
Shipping	Ice bag