


Anti-Parp2 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Ugn5

Item Number: CM0019772



Introduction

Rabbit polyclonal antibody targeting human PARP2, suitable for Western blot and ELISA applications. Cross-reacts with mouse and rat. Ideal for DNA repair, chromatin remodeling, and ADP-ribosylation studies. 66kDa.

[Learn More](#)

Parameter	Value
Product Name	PARP2 Rabbit pAb
Remarks/Alias	ARTD2; ADPRT2; PARP-2; ADPRTL2; ADPRTL3; pADPRT-2; PARP2
Species	Human
Gene ID (Human)	10038
Gene ID	10038
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9UGN5
Protein Weight	66kDa
Shipping	Ice bag