

Anti-Parp4 Polyclonal Antibody For Wb, Elisa - Q9Ukk3

Item Number: CM0026860



Introduction

Rabbit polyclonal antibody against human PARP4, validated for WB and ELISA. Detects PARP4 in human, mouse, and rat samples. Ideal for vault complex and ADP-ribosylation research.

[Learn More](#)

Parameter	Value
Product Name	PARP4 Rabbit pAb
Remarks/Alias	PH5P; p193; ARTD4; PARPL; VPARP; VWASC; PARP-4; VAULT3; ADPRTL1; PARP4
Species	Human
GeneID (Human)	143
GeneID	143
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human PARP4 (NP_006428.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UKK3
Protein Weight	193kDa
Shipping	Ice bag