

Anti-Parp4 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Ukk3

Item Number: CM0022485



Introduction

Rabbit polyclonal antibody against human PARP4, validated for Western blot and ELISA. Reacts with human, mouse, and rat. Ideal for vault particle and ADP-ribosylation studies.

[Learn More](#)

Parameter	Value
Product Name	PARP4 Rabbit pAb
Remarks/Alias	PH5P; p193; ARTD4; PARPL; VPARP; VWA5C; PARP-4; VAULT3; ADPRTL1; PARP4
Species	Human
GeneID (human)	143
GeneID	143
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UKK3
Protein Weight	193kDa
Shipping	Ice bag