

Anti-Parp8 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8N3A8

Item Number: CM0018669

Introduction

PARP8 rabbit polyclonal antibody validated for WB, IHC-P, ELISA in Human, Mouse, Rat. Targets protein mono-ADP-ribosyltransferase PARP8 (ARTD16), involved in ADP-ribosylation. Ideal for cancer and cell signaling research.

[Learn More](#)

Parameter	Value
Product Name	PARP8 Rabbit pAb
Remarks/Alias	ARTD16; pART16; PARP8
Species	Human
Gene ID (Human)	79668
Gene ID	79668
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human PARP8 (NP_001171527.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N3A8
Protein Weight	96kDa
Shipping	Ice bag