

Anti-Parva Polyclonal Antibody For Wb And Elisa - Q9Nvd7

Item Number: CM0018738



Introduction

Rabbit polyclonal antibody against human PARVA (alpha-parvin), validated for Western blot and ELISA. Cross-reacts with mouse and rat. Suitable for cell adhesion, cytoskeleton, and cardiovascular research.

[Learn More](#)

Parameter	Value
Product Name	PARVA Rabbit pAb
Remarks/Alias	PARVA; CH-ILKBP; MXRA2; alpha-parvin
Species	human
Gene ID	55742
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human PARVA (NP_060692.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NVD7
Protein Weight	20kDa/42kDa
Shipping	Ice bag