

Anti-Parvalbumin (Pvalb) Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P20472

Item Number: CM0016753



Introduction

Rabbit polyclonal antibody targeting human Parvalbumin alpha (PVALB). Validated for WB, IF-P, IHC-P, ELISA; cross-reacts with human, mouse, rat. Ideal for muscle physiology and calcium-binding protein research.

[Learn More](#)

Parameter	Value
Product Name	Parvalbumin/PVALB Rabbit pAb
Remarks/Alias	D22S749; Parvalbumin (PVALB)
Species	Human
GeneID (Human)	5816
GeneID	5816
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-110 of human Parvalbumin (NP_002845.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20472
Protein Weight	12kDa
Shipping	Ice bag