

Anti-Parvalbumin/Pvalb Polyclonal Antibody For Wb - P20472

Item Number: CM0025637



Introduction

Rabbit polyclonal antibody targeting Parvalbumin/PVALB (alpha-PV). Suitable for WB, IF-P, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Used to investigate calcium-dependent muscle relaxation and calcium signaling.

[Learn More](#)

Parameter	Value
Product Name	Parvalbumin/PVALB Rabbit pAb
Remarks/Alias	D22S749; Parvalbumin (PVALB)
Species	Human
Gene ID (Human)	5816
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-110 of human Parvalbumin (PVALB) (NP_002845.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20472
Protein Weight	12kDa
Shipping	Ice bag