

Anti-Parvalbumin/Pvalb Rabbit Monoclonal Antibody For Wb, Ihc-P, If-P, Elisa - P20472

Item Number: CM0002764



Introduction

Rabbit monoclonal anti-Parvalbumin/PVALB antibody validated for WB, IF-P, IHC-P, ELISA. Suitable for human, mouse, rat samples. Recognizes a ~12 kDa calcium-binding protein involved in muscle relaxation.

[Learn More](#)

Parameter	Value
Product Name	Parvalbumin/PVALB Rabbit mAb
Remarks/Alias	D22S749; Parvalbumin (PVALB)
Species	Human
GeneID (human)	5816
GeneID	5816
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human Parvalbumin (PVALB) (P20472)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20472
Protein Weight	12kDa
Shipping	Ice bag