


Anti-Parvb Rabbit Polyclonal Antibody For Wb And Elisa - Q9Hbi1

Item Number: CM0018251



Introduction

Rabbit polyclonal antibody targeting human PARVB (beta-parvin, affixin). Validated for Western blotting and ELISA, with cross-reactivity to mouse and rat. Ideal for studying cell adhesion, migration, and integrin signaling in muscle and other tissues.

[Learn More](#)

Parameter	Value
Product Name	PARVB Rabbit pAb
Remarks/Alias	CGI-56; PARVB
Species	Human
GeneID	29780
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human PARVB (NP_037459.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HBI1
Protein Weight	42kDa
Shipping	Ice bag