


Anti-Pask Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q96Rg2

Item Number: CM0026395



Introduction

Recombinant-based rabbit polyclonal antibody targeting human PASK kinase (STK37/PASKIN). Validated for WB, IHC-P, IF/ICC, ELISA with human, mouse, rat reactivity. Ideal for energy homeostasis and translational signaling research.

[Learn More](#)

Parameter	Value
Product Name	PASK Rabbit pAb
Remarks/Alias	STK37; PASKIN; PASK
Species	Human
GeneID	23178
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1194-1323 of human PASK (NP_055963.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96RG2
Protein Weight	143kDa
Shipping	Ice bag