

Anti-Patz1 Monoclonal Antibody For Wb, Elisa - Q9Hbe1

Item Number: CM0011023



Introduction

Rabbit monoclonal antibody targeting human PATZ1 (POZ-, AT hook-, and zinc finger-containing protein 1), a transcriptional regulator. Validated for Western blot and ELISA. Suitable for studies in T-cell development, senescence, and cancer.

[Learn More](#)

Parameter	Value
Product Name	PATZ1 Rabbit mAb
Remarks/Alias	ZSG; MAZR; PATZ; RIAZ; ZBTB19; ZNF278; dj400N23
Species	Human
Gene ID	23598
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 250-350 of human PATZ1 (NP_055138.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9HBE1
Protein Weight	74kDa
Shipping	Ice bag