

Anti-Patz1 Polyclonal Antibody For Wb, Elisa - Q9Hbe1

Item Number: CM0018083



Introduction

High-quality rabbit polyclonal antibody against human PATZ1 for Western blot and ELISA applications. Specificity for human, mouse, and rat samples. Suitable for transcription and cancer research.

[Learn More](#)

Parameter	Value
Product Name	PATZ1 Rabbit pAb
Remarks/Alias	ZSG; MAZR; PATZ; RIAZ; ZBTB19; ZNF278; dj400N23; PATZ1
Species	Human
GeneID (Human)	23598
GeneID	23598
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 250-350 of human PATZ1 (NP_055138.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HBE1
Protein Weight	74kDa
Shipping	Ice bag