

Anti-Pax1 Polyclonal Antibody For Wb And Elisa - P15863

Item Number: CM0024461

Introduction

Rabbit polyclonal antibody targeting human PAX1, a developmental transcription factor. Suitable for Western blot and ELISA. Ideal for embryogenesis and vertebral column research.

[Learn More](#)

Parameter	Specification
Product Name	PAX1 Rabbit pAb
Remarks/Alias	OFC2; HUP48; PAX1
Species	Human
Gene ID (Human)	5075
Gene ID	5075
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 300-400 of human PAX1 (NP_006183.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P15863
Protein Weight	55kDa
Shipping	Ice bag