

Anti-Pax2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q02962

Item Number: CM0029282

Introduction

Rabbit polyclonal antibody targeting human PAX2, validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for kidney development and urogenital tract research.

[Learn More](#)

Parameter	Value
Product Name	PAX2 Rabbit pAb
Remarks/Alias	FSGS7; PAPRS; PAX-2; PAX2
Species	Human
GeneID (Human)	5076
GeneID	5076
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 150-320 of human PAX2 (NP_003978.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02962
Protein Weight	45kDa
Shipping	Ice bag