

Anti-Pax3 Rabbit Polyclonal Antibody For Wb And Elisa - P23760

Item Number: CM0016108



Introduction

High-quality rabbit polyclonal antibody against human PAX3, validated for western blotting (WB) and ELISA. Reacts with human, mouse, and rat samples. Ideal for research on neural development, myogenesis, and transcription regulation.

[Learn More](#)

Parameter	Value
Product Name	PAX3 Rabbit pAb
Remarks/Alias	WS1; WS3; CDHS; HUP2; PAX-3; PAX3
Species	Human
GeneID (Human)	5077
GeneID	5077
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23760
Protein Weight	53kDa
Shipping	Ice bag