

Anti-Pax3 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P23760

Item Number: CM0019638



Introduction

Rabbit polyclonal antibody against PAX3, validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Ideal for research on neural development, myogenesis, and Waardenburg syndrome.

[Learn More](#)

Parameter	Value
Product Name	PAX3 Rabbit pAb
Remarks/Alias	WS1; WS3; CDHS; HUP2; PAX-3; PAX3
Species	Human
Gene ID (Human)	5077
Gene ID	5077
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-236 of human PAX3 (NP_852122.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23760
Protein Weight	53 kDa
Shipping	Ice bag