

Anti-Pax5 Rabbit Monoclonal Antibody For Wb, Ihc-P, Ip, And Elisa - Q02548

Item Number: CM0004971



Introduction

Rabbit monoclonal antibody targeting human PAX5 (BSAP), suitable for WB, IHC-P, IP, and ELISA. Reliable tool for B-cell lineage commitment and lymphoma studies.

[Learn More](#)

Parameter	Value
Product Name	PAX5 Rabbit mAb
Remarks/Alias	ALL3; BSAP; PAX-5; PAX5
Species	Human
Gene ID (Human)	5079
Gene ID	5079
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human PAX5 (Q02548).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human
SWISS	Q02548
Protein Weight	42kDa
Shipping	Ice bag