

Anti-Pax6 Monoclonal Antibody For Wb/Ihc-P/Elisa - P26367

Item Number: CM0003345



Introduction

High-quality PAX6 rabbit monoclonal antibody validated for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat PAX6, a key transcription factor in eye and CNS development. Ideal for IVD and research applications.

[Learn More](#)

Parameter	Value
Product Name	PAX6 Rabbit mAb
Remarks/Alias	AN; PAX6
Species	Human
GeneID (Human)	5080
GeneID	5080
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 323-422 of human PAX6 (NP_000271.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P26367
Protein Weight	47kDa/48kDa
Shipping	Ice bag