

Anti-Pax6 Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - P26367

Item Number: CM0030800



Introduction

PAX6 rabbit polyclonal antibody validated for WB, IF/ICC, IP, ELISA in human, mouse, rat. Targets the paired box transcription factor critical for eye, CNS, and pancreas development. Ideal for developmental biology and disease research.

[Learn More](#)

Parameter	Value
Product Name	PAX6 Rabbit pAb
Remarks/Alias	AN; AN1; AN2; FVH1; MGDA; WAGR; ASGD5; D11S812E; PAX6
Species	Human
Gene ID (Human)	5080
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human PAX6 (NP_000271.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P26367
Protein Weight	47kDa
Shipping	Ice bag