

Anti-Pax6 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P26367

Item Number: CM0007503



Introduction

Rabbit monoclonal antibody targeting PAX6, validated for WB, IHC-P, IF/ICC, and ELISA. Recognizes human, mouse, and rat PAX6. Ideal for developmental biology and aniridia research.

[Learn More](#)

Parameter	Value
Product Name	PAX6 Rabbit mAb
Remarks/Alias	AN; PAX6
Species	Human
GeneID	5080
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 323-422 of human PAX6 (NP_000271.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P26367
Protein Weight	47kDa/48kDa
Shipping	Ice bag