

Anti-Pax7 Rabbit Polyclonal Antibody For Wb, Elisa - P23759

Item Number: CM0024739

Introduction

Rabbit polyclonal antibody targeting human PAX7 (HuP1). Validated for WB and ELISA. Cross-reacts with mouse and rat. Suitable for muscle stem cell and myogenesis research. ~55 kDa.

[Learn More](#)

Parameter	Value
Product Name	PAX7 Rabbit pAb
Remarks/Alias	HUP1; RMS2; PAX7B; CMYP19; MYOSCO; PAX7
Species	Human
GeneID (Human)	5081
GeneID	5081
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 350-470 of human PAX7 (NP_039236.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23759
Protein Weight	55kDa
Shipping	Ice bag