


Anti-Pax8 Rabbit Monoclonal Antibody For If/Icc, Ihc-P, Elisa - Q06710

Item Number: CM0007530



Introduction

High-quality rabbit monoclonal antibody targeting human PAX8, validated for IF/ICC, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for thyroid transcription factor studies and Wilms tumor research.

[Learn More](#)

Parameter	Value
Product Name	PAX8 Rabbit mAb
Remarks/Alias	PAX-8; PAX8
Species	Human
GeneID	7849
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human PAX8 (NP_003457.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q06710
Protein Weight	48kDa
Shipping	Ice bag