


Anti-Pax8 Rabbit Polyclonal Antibody For Wb & Elisa - Q06710

Item Number: CM0014272



Introduction

High-quality PAX8 Rabbit Polyclonal Antibody validated for WB and ELISA. Suitable for human, mouse, and rat samples. Ideal for thyroid transcription factor research and developmental biology studies. Backed by UniProt Q06710.

[Learn More](#)

Parameter	Value
Product Name	PAX8 Rabbit pAb
Remarks/Alias	PAX-8; PAX8
Species	Human
GeneID	7849
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 200-300 of human PAX8 (NP_003457.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q06710
Protein Weight	48kDa
Shipping	Ice bag